
Ship P. Counter

Ship P. Counter enables to obtain the number of phytoplankton continuously by automatic inspection monitoring system and receive the result from a sample of ballast water

Features

- Automation & Reduced inspection time
- Cost reduction
- Fast: Results are delivered in 3 minutes
- High sensitivity: Counting Individual Living Cell
- No dye-substrates
- Easy sampling & counting



* Sample image


IMO Agenda Proposal



INTERNATIONAL MARITIME ORGANIZATION

MEPC 74-INF.18 agenda proposal





Ballast water, representative medium of exotic species inflow

Fully connected Convolutional Convolutional laver layer 1 layer 2 Max pooling laver 2 Output Max pooling laver 1 lavers Input laye

Deep Learning Algorithm for Inspection

Detecting algorism




1.Images Filter Process

2. Deep Learning Algorithm Process





Specification



Ship P.Counter		
Interrogated Volume	1.4mL * 8times	
Sensitivity	< 1cell / mL	
Dynamic Range	0-1000cell / mL	
Time to result	< 3minutes for 1Cycle	
Power	AC 220V / 50, 60Hz	
Connectivity	USB or Ethernet or LTE(option)	
Dimensions	700 * 700 * 400 mm	
Operating Temperature	0 ~ 50°C	
Storage Temperature	-10 ~ 55°C	

* The above specifications are subject to change depending on the company's circumstances

Sequence of operation

- 1. Drain water from inside
- 2. Open valve 1 to store water in the water tank
- 3. Open valves 3 and 4 for air discharge at the rear end of the water tank when the water tank reaches the intermediate level (from 1 to 1.5L)
- 4. Close all valves and start filming
- 5. Drain the area (600ml~1.2L) and close the valve again.
- 6. Repeat movement 4 and 5
- 7. Drain all water after 8 shots
- 8. Go back to the movement 2 and repeat the whole thing

Switch & Mode

Switch Status		Mode
OFFON'S	Valve Switch OffControl Unit Switch Off	Power Off
XOFFERING	Valve Switch ONControl Unit Switch Off	Water Circulation (Cleaning) Mode
	Valve Switch ONControl Unit Switch ON	Automatic Check Mode

